

Fluorescence *in situ* Hybridization Protocol 60°C for use with KromaTiD Pinpoint FISH™ Probes

Reagents & Materials Provided

- FISH Hybridization Buffer (provided by KromaTiD)
- Probe Solution (provided by KromaTiD)

Samples

- This protocol requires use of cell samples that have been fixed in 3:1 MeOH: Glacial acetic acid (Carnoy's fixative), prepared on cytogenetic-quality microscope slides (See Appendix A).

Reagents Required

- 20X Saline Sodium Citrate (SSC) Buffer (pH7.0) (see Appendix A)
- 2X SSC Buffer (pH7.0) (see Appendix A)
- 0.4X SSC / 0.05% Tween-20 (see Appendix A)
- Antifade mounting medium with DAPI

Equipment and Supplies

- Glass slides (25x75x1mm, positively charged)
- Coverslips (22x50mm)
- 50mL Coplin jars (acetal or glass)
- Water bath (60°C)
- Slide hybridization system (Cytobrite or equivalent hot plate)
- Incubator (37°C)
- Humidified chamber (see Appendix A)
- Calibrated thermometer
- Rubber Cement
- Vortex mixer
- Benchtop microcentrifuge
- Micropipette and sterile tips (2-200µL range)
- Microcentrifuge tubes (sterile, DNase/RNase free)
- Forceps

Fluorophore Information

Fluorophore	Excitation/Emission (nm)
Atto425 / Aqua	436 / 485
6-FAM / Spectrum Green	490 / 525
Atto550 / Spectrum Orange	555 / 576
TexRed	595 / 620
Atto643 / 647 / Cy5	643 / 669

Handling/Storage Information

- Store all probes and hybridization buffers at -20°C upon receiving.
- Probes are to be protected from light until ready to begin.

Pinpoint FISH™ Protocol

Slide Preparation for Pinpoint FISH™

1. Samples for Pinpoint FISH™ hybridization should be prepared using standard metaphase spread preparation techniques on glass slides specified for cytogenetic applications (positively charged).
 - a. Submerging new slides in 100% ethanol may help remove any residue that could be present from the slide manufacturing process.
 - b. For good metaphase spreads, drying should take place in a humidity-controlled chamber to optimize chromosome spreading.
 - c. Freshly dropped slides may be artificially aged for 10 minutes in a 90°C oven.

Hybridization Procedure

Note: Prepare the hybridization humidity chamber with DI H₂O and place in the 37°C incubator to equilibrate prior to beginning.

2. Prepare the hybridization mixture for one full slide (30µL), or in a 3:2 ration of FISH Hybridization Buffer to probe solution if a smaller volume of hybridization mix is desired.

Note: Keep probe products on ice and away from full light if not in use, preferably in the dark.

Solution	Volume
FISH Hybridization Buffer	18 µL
Probe Solution	12 µL

Note: If combining more than one KromaTiD Pinpoint FISH probe in a single assay, then combine 12uL of each probe into a single tube and proceed to dehydrate the probes. Resuspend the dried probe mixture in 12uL of nuclease-free H₂O.

- a. Vortex the mixture three times for 5 seconds each.
 - b. Centrifuge for 10 seconds.
3. Apply hybridization mix to the slide using a micropipette. Use care not to introduce air bubbles.
4. Apply a coverslip
 - a. If a bubble is present underneath the coverslip, it can be forced out with gentle and even pressure to the coverslip.

Note: Pressing on a coverslip can physically damage cells if not done with great care.

5. Seal the edges of the coverslip with rubber cement. Be sure all edges are sealed to prevent the hybridization mixture from drying out.

Note: Ensure that no rubber cement gets underneath the slide. It can obstruct slide contact with the heat block in the denature step that follows.

6. Denature the slide using a hybridization system (Cytobrite or equivalent hot plate) by placing the slide directly onto the heat source at 75°C for 5 minutes.
7. After the program ends, transfer the slide to a pre-warmed humidified chamber. Place the humidified chamber containing the hybridization slides into the 37°C incubator for a minimum of 4 hours.

Note: An overnight hybridization can be performed in a humidified chamber for up to 20 hours at 37°C.

Post-Hybridization Slide Treatment

1. Prepare wash solution 1: place 50mL of 0.4X SSC/ 0.05% Tween-20 solution in a Coplin jar and pre-heat in a water bath set to 60°C.
2. Prepare wash solution 2: place 50mL of 2X SSC solution in a Coplin jar and leave at room temperature.
3. Gently roll off the dry rubber cement seal and remove coverslip. The coverslip can be removed by submersion in a 2X SSC Coplin jar at room temperature, or the coverslip can be removed gently by hand.
4. Ensure the temperature of the 0.4X SSC/ 0.05% Tween-20 (wash solution 1) inside the Coplin is at 60°C. Add the slides and wait 1 minute.
5. Move the slides to wash solution 2 (2X SSC at RT).
6. After 1-2 minutes, remove the slides from wash solution 2 and airdry in the dark.
***Note:** Wipe the excess solution from the back of the slide using a lint-free towel or equivalent will help to reduce background during imaging. Also, propping slides up on their sides or label-end down will prevent residues such as ink or adhesive potentially running down onto the sample from the label area.*
7. Add 18µL of Antifade mounting medium with DAPI. Apply a coverslip with care to avoid trapping air bubbles.
8. Allow the DAPI stain to react for at least 15 minutes prior to imaging.

Appendix A: Reagents

20x Saline Sodium Citrate (SSC)

Reagent	Amount
Sodium Chloride	175.3g
Sodium Citrate	88.2g
DI H ₂ O	QS to 1L
pH to 7.0	

2x Saline Sodium Citrate (SSC)

Reagent	Amount
20x SSC	100mL
DI H ₂ O	QS to 1L
pH to 7.0	

0.4x SSC / 0.05% Tween-20

Reagent	Amount
20X SSC (pH 7.0)	20mL
Tween-20	500μL
DI H ₂ O	QS to 1L

Cytogenetic-Quality Microscope Slides

SciGene CytoDrop® Slides Cat# 2060-01-01 or equivalent

Humidified Chamber

KromaTiD uses 96-well plastic freezer box with the outside wells filled with DI H₂O as the humidified chamber. Slides are placed flat and upright in the chamber before sealing the box with parafilm and placing the entire chamber into a humidified-controlled 37°C incubator.

Customer Notification

1. All Products and Deliverables are supplied for internal scientific research purposes only and are not intended for i) human consumption, including, but not limited to, foods or pharmaceuticals, ii) diagnostic purpose including, but not limited to, human or veterinary *in vivo* or *in vitro* diagnostics, or use in cosmetics or other goods. Research purposes means *in vitro* laboratory studies or *in vivo* use in laboratory organisms only.
2. Products and Deliverables are supplied for the Customer's personal research activities and noncommercial use. Products and Deliverables must not be sold or otherwise redistributed without KromaTiD's consent.
3. KromaTiD's full Terms of Sale can be found on our website.



TECHNICAL SUPPORT:
LIVE CHAT: M-F 9 AM- 5 PM MT
EMAIL: TECHSUPPORT@KROMATID.COM

INQUIRIES AND ORDERS:
EMAIL: SALES@KROMATID.COM
WEB: KROMATID.COM
Ph: 720.815.2898
Fx: 720.815.2902

For Research Use Only. Not for use in diagnostic procedures.